

WVSP TOXICOLOGY BLOOD DRUG PANEL

A negative reporting statement can signify that no compounds were detected in the specimen, reporting criteria was not met or that the Limit of Detection/Limit of Quantitation reporting threshold was not met. Levels listed on this Panel are informational and should not be used for interpretation of quantitative results listed on the report.

	P1-Panel 1 P2-Panel 2 N2- Negative 2 Panel T-THC/Metabolite Panel	Compound	Screen Limit of Detection (ng/ml)	Limit of Quantitation (ng/ml)
1	P1	2-hydroxyethylflurazepam	5	5
2	P1	6-MAM	0.5	0.5
3	P1	7-aminoclonazepam	5	5
4	P1	7-aminonitrazepam	5	5
5	P1	Acetyl Fentanyl	0.5	0.5
6	P1	Acetyl Norfentanyl	0.5	0.5
7	P1	Alprazolam	5	5
8	P1	Amphetamine	5	5
9	P1	Benzoyllecgonine	5	5
10	P1	Buprenorphine	0.5	1
11	P1	Carisoprodol	200	200
12	P1	Clonazepam	5	5
13	P1	Cocaine	5	5
14	P1	Codeine	5	5
15	P1	Cyclobenzaprine	5	5
16	P1	Dextromethorphan	5	5
17	P1	Dextrorphan	5	5
18	P1	Diazepam	5	5
19	P1	Dihydrocodeine	5	5
20	P1	EDDP	5	5
21	P1	Fentanyl	0.5	0.5
22	P1	Flunitrazepam	5	5
23	P1	Flurazepam	5	5
24	P1	Hydrocodone	5	5
25	P1	Hydromorphone	5	5
26	P1	Lorazepam	5	5
27	P1	3,4-Methylenedioxyamphetamine (MDA)	5	5
28	P1	3,4-Methylenedioxy-N-ethylamphetamine (MDEA)	5	5
29	P1	3,4-Methylenedioxymethamphetamine (MDMA)	5	5
30	P1	Meprobamate	200	200
31	P1	Methadone	5	5
32	P1	Methamphetamine	5	5
33	P1	Midazolam	5	5
34	P1	Morphine	5	5
35	P1	Nitrazepam	5	5
36	P1	Norbuprenorphine	0.5	0.5
37	P1	Nordiazepam	5	5
38	P1	Norfentanyl	0.5	0.5
39	P1	O-Desmethyltramadol	5	5
40	P1	Oxazepam	10	10
41	P1	Oxycodone	5	5
42	P1	Oxymorphone	5	5
43	P1	Phenazepam	5	5
44	P1	Phencyclidine	5	5

WVSP TOXICOLOGY BLOOD DRUG PANEL

	P1-Panel 1 P2-Panel 2 N2- Negative 2 Panel T-THC/Metabolite Panel	Compound	Screen Limit of Detection (ng/ml)	Limit of Quantitation (ng/ml)
45	P1	Phentermine	5	5
46	P1	Sufentanil	0.5	0.5
47	P1	Tapentadol	5	5
48	P1	Temazepam	5	5
49	P1	Tramadol	5	5
50	P1	Zaleplon	5	5
51	P1	Zolpidem	5	5
52	P1	Zopiclone	5	5
53	P1	α -hydroxyalprazolam	5	5
54	P2	7-aminoflunitrazepam	5	5
55	P2	Bupropion	5	5
56	P2	Carbamazepine	5	5
57	P2	Chlordiazepoxide	5	5
58	P2	Citalopram	5	5
59	P2	Desalkylflurazepam	5	5
60	P2	Desmethyldoxepin	5	5
61	P2	Desmethyilmirtazapine	5	5
62	P2	Diphenhydramine	5	5
63	P2	Doxepin	5	5
64	P2	Doxylamine	5	5
65	P2	Estazolam	5	5
66	P2	Etizolam	5	5
67	P2	Fluoxetine	5	5
68	P2	Ketamine	5	5
69	P2	Meperidine	5	5
70	P2	Methylphenidate	5	5
71	P2	Mirtazapine	5	5
72	P2	Norephedrine	5	5
73	P2	Norfluoxetine	-	5
74	P2	Norketamine	5	5
75	P2	Normeperidine	5	5
76	P2	Pentazocine	5	5
77	P2	Promethazine	5	5
78	P2	Pseudoephedrine/Ephedrine	5	5
79	P2	Quetiapine	5	5
80	P2	Sertraline	5	5
81	P2	Trazodone	5	5
82	P2	Venlafaxine	5	5
83	N2	Butabarbital	200	200
84	N2	Butalbital	200	200
85	N2	Pentobarbital/Amobarbital	200	200
86	N2	Phenobarbital	200	200
87	N2	Secobarbital	200	200
88	T	Delta-9 Tetrahydrocannabinol (THC)	1	1
89	T	11-Nor-9-carboxy- Δ -9 tetrahydrocannabinol (THC-COOH)	5	5
90	T	11-hydroxy- Δ -9 tetrahydrocannabinol (THC-OH)	1	1